

# Specifications TV12LT

Tamson Visibility bath 15 liters



- ⊕ Small footprint
- ⊕ Detachable front window
- ⊕ Internal LED light
- ⊕ Ultra high stability
- ⊕ RS232 communication
- ⊕ Bath drain
- ⊕ Built-in cooling
- ⊕ -50°C / -40°F
- ⊕ 4 places, small bath volume

| Item              | Unit     | TV12LT                                |             |            |
|-------------------|----------|---------------------------------------|-------------|------------|
| <b>Range*</b>     |          | -40..+20°C                            | -40..+80°C  | -50..+20°C |
| <b>ordering</b>   |          | -40..+68°F                            | -40..+176°F | -58..+68°F |
| 230V/50Hz         |          | 00T0410                               | 00T0425     | 00T0470    |
| 115V/60Hz         |          | 00T0415                               | 00T0430     | 00T0475    |
| 230V/60Hz         |          | 00T0420                               | 00T0435     | 00T0480    |
| Reading           | °C or °F | menu selectable                       |             |            |
| Interface         |          | RS232                                 |             |            |
| Setting ±         | [°C/F]   | 0.01                                  |             |            |
| Stability **      | [°K]     | stdev 0.004<br>min/max ±0.014         |             |            |
| Uniformity **     | [°K]     | stdev 0.004<br>min/max ±0.013         |             |            |
| Heating           | [W]      | 500 + 700                             |             |            |
| Heaters           |          | 2                                     |             |            |
| Bath volume       | [L]      | 15                                    |             |            |
| Number of lids    |          | 1 lid with 4 x round diameter<br>51mm |             |            |
| Window            | [mm]     | 255x230                               |             |            |
| Opening bath      | [mm]     | 248*73                                |             |            |
| Depth             | [mm]     | 300                                   |             |            |
| Length            | [mm]     | 670                                   |             |            |
| Width             | [mm]     | 425                                   |             |            |
| Height            | [mm]     | 720                                   |             |            |
| Weight            | [kg]     | 65                                    |             |            |
| Power             | [Watt]   | Nominal 800<br>Maximum 2100           |             |            |
| Ambient condition | [°C]     | 18 .. 23                              |             |            |
| CE                |          | All models conform to CE regulation   |             |            |

\* +/- 2°C production variation  
 \*\* Measured @-20°C in methanol

## General

Tamson viscometer and Tamson calibration baths are specially designed for tests that require ultra precise temperature control, or processes that need to be followed visually, e.g. viscometry (conforms to ASTM D445, IP71-1), thermometer and sensor calibration, density and reaction rate measurement, etc. The bath is fitted with a double window of which the front pane is detachable for cleaning purposes. The window is heated preventing built up of condensate.

## Construction

The stainless steel construction ensures exceptionally stable temperatures which is further improved by an ingenious stirring mechanism with baffle plates. All wetted parts are made of stainless steel and PTFE, providing resistance against all usual bath fluids. The bath is fitted with adjustable feet for levelling. The cover of the bath has 4 round 51 mm holes with lids, for suspending glass capillary viscometers in holders. The use of thermo insulated windows and window heating ensures clean sight. A permanent LED light is located in the top plate to supply clear light and guarantee optimal visibility inside the bath. A bath overflow outlet protects against expanding bath oil when the bath filling is too high.

## Agitation

A vane type stirrer with maintenance free bearings moves the bath fluid past a special heater ensuring optimal control temperature and excellent uniformity.

## Span

Dependable from ordering code. Span lies from -50°C/-58°F up to + 80°/176°F

## Safety

The bath conforms to CE-regulation. Further the bath is equipped with a mechanical over temperature device which trips when in case of malfunction the bath exceeds the preset maximum temperature. This feature guarantees safe around the clock operation.

## Fine adjustment and offset

After the bath has stabilized the set point may be more accurately adjusted in the range of -5.00° to + 5.00°, if necessary.

## Ambient condition

For proper cooling performance ambient temperature needs to be within the range of 18 .. 23°C.

## Options

- Optical Level indicator **07T0080**

## Accuracy

Recovery from temperature dip


TV12 bath



Conventional bath



Inside glass viscometer tube




## Accuracy

Methanol

standard deviation  $\pm 0.004^\circ\text{C}$

min / max  $\pm 0.013^\circ\text{C}$



## Homogeneity


Methanol

standard deviation  $\pm 0.004^\circ\text{C}$


min / max  $\pm 0.013^\circ\text{C}$

- 1: Stirrer
- 2: Heater
- 3: Baffle plate
- 4: Circulation
- 5: Cooling















## TV12LT




## Conventional bath



## Temperature gradient TV12LT versus conventional system

| Item  |   |   |              |                  |
|---|---|---|--------------|------------------|
| Timer   |    | 10T6090   |              |                  |
| Bath fluid  |    | See datasheet "Bath fluids"   |              |                  |
| Thermometers  |    | ASTM nr.  | Ordering no. | Range°C          |
|   |   | 44C   | 25T0937      | +18.5 ... +21.5  |
|   |   | 46C   | 25T0938      | +48,6 .. +51,4   |
|   |   | 120C  | 25T0990      | +38.6 ... +41.4  |
|   |   | 46C   | 25T0939      | +48.6 ... +51.4  |
|   |   | 47C   | 25T0940      | + 58.6 ... +61.5 |
|   |   | 121C  | 25T0991      | +98.6 ... +101.4 |
|   |   | 127C  | 25T0995      | -21,4 ... -18,6  |
|   |   | 128C  | 25T0994      | -1,4 .. 1,4      |
|   |   | Other ranges available on request   |              |                  |
| Thermometer holder  |   | 00T0239   |              |                  |
| Calibration,reference oils                                      |    | See datasheet "Viscosity calibration standards"                                     |              |                  |
| Glass viscometers   |   | See datasheet "Viscometer to ASTM D446, IP71 and BS188"                             |              |                  |
| Viscometer holders  |  | See datasheet "Viscometer holders"  |              |                  |
| Cover lid   |  | This item is standard included  | 13T3006.08   |                  |
| Cover lid   |  | Cover with 4 holes of 60 mm for Pinkevitch viscometers                              | 03T2117      |                  |
| Cover lid   |  | Optional to be ordered for 4 full position Cannon Fenske Transparent                | 13T3006.09   |                  |
| Cover lid, improved contrast viscometer reading with 13T3006.08 |  |  | 13T6220      |                  |
| Leveling platform   |  |  | 13T6200      |                  |
| Calibration/metall block  |  |   | 13T6210      |                  |


## Levelling platform(13T6220) with high background contrast




4 opaline glass plates are mounted under the top lid.

The semi transparent white glass realises uniform background and optimizes contrast and readout of the viscometer.

## Levelling platform(13T6200) and metal block(13T6210)



## Dimensions top plate



☞ Number of holes, dimensions and position  
Other dimensions on request ☞