


# CCTO A50T950

Non thermostated  
Effective volume: 150 ml to 300 ml  
For large flat sample  
Sample active  $\varnothing$ : 7.5 cm  
Sample active area: 45 cm<sup>2</sup>



# CCTO

## Large Flat Sample Cell

### A50T950

#### Description

The CCTO is a pyrex glass body cell designed for testing corrosion on large flat samples. The shape of the cell combined to the large area of the working electrode and the stainless steel counter electrode make it particularly suitable for testing low conductivity media and/or highly resistive samples.

#### Ordering Information

The CCTO product includes

<i>Part No.</i>	<i>Qty</i>	<i>Type</i>	<i>Description</i>
A50T966	1	CCTO/EM	CCTO cell bottom
A50T963	1	CCTO/CC	Pyrex glass body for CCTO cells
A50T967	1	CCTO/CV	CCTO cell cap
A50T968	1	CCTO/EA	Stainless steel auxiliary electrode for CCTO cells
A50T971	1	CCTO/V15/J8	Ground joint for auxiliary electrode (CCTO)
A50T965	1	CCTO/PJ	Salt-bridge extension for CCTO cells
A50T975	1	CCTO/TG	Gas inlet tubing for CCTO cells
A50T970	1	CCTO/TR	Gas outlet tubing for CCTO cells

#### CCTO/TH consumables

<i>Part No.</i>	<i>Qty</i>	<i>Type</i>	<i>Description</i>
X31V001	1		Silicone washer for CCTO cell (4 pcs)

D11V010 © Radiometer Analytical SAS • France • 2007-02D. All rights reserved. Data subject to change without notice



**RADIOMETER ANALYTICAL SAS**

72 rue d'Alsace, 69627 Villeurbanne Cedex, France

E-mail: [radiometer@analytical.com](mailto:radiometer@analytical.com) Web: [www.radiometer-analytical.com](http://www.radiometer-analytical.com)

Tel. : +33 (0)4 78 03 38 38 - Fax : +33 (0)4 78 68 88 12