

Mantas calefactoras | Heating mantles

Marca|Trademark: Nahita

Las mantas calefactoras con regulación electrónica de temperatura son unos instrumentos básicos en cualquier laboratorio para procesos de destilación, digestión, evaporación, extracción y ebullición. Con mueble exterior de acero pintado y aislamiento térmico interior por lana de fibra mineral, disponen de elementos calefactores de cromo-níquel. De calefacción rápida y homogénea y consumo reducido, incluyen un indicador luminoso de conexión y lámpara de señalización de funcionamiento. Temperaturas hasta 370 °C cuentan con un dispositivo para sujetar varillas.

Heating mantles with electronic temperature regulation are essential equipments in any laboratory for distillation, digestion, evaporation, extraction and boiling. The outer frame is made of furnace-painted steel. The inner vessel is thermally isolated with mineral fiber and the heating elements are composed of chromium-nickel. Quick and homogeneous heating and low consumption, they include a connection light device and on/off lamp. For temperatures up to 370 °C. Clamp for supporting rods included.



Referencia	Capacidad	Diámetro int.	Consumo	Dimensiones
Code	Capacity	Inner diameter	Consumption	Size
50655210	100 mL	70 mm	100 W	220x160 mm
50655220	250 mL	90 mm	150 W	160x195 mm
50655250	500 mL	105 mm	250 W	220x160 mm
50655310	1000 mL	150 mm	350 W	300x210 mm
50655320	2000 mL	175 mm	450 W	220x225 mm
50655330	3000 mL	188 mm	600 W	320x230 mm
50655350	5000 mL	240 mm	800 W	325x245 mm