

# ANDREESCU LABOR & SOFT SRL

Stoc produse din sticla

In prezenta lista veti gasi urmatoarele produse : Cilindri, baloane cotate, pipete gradate, Exsicatoare, pipete cotate, biurete, pahare Berzelius gradate, pahare Erlenmeyer, pahare Erlenmeyer pentru vid, palnii, baloane Kjeldhal, baloane cu fund plat, sticle brune cu dop rodat, eprubete gradate, sticle de ceas ...



Preturile sunt exprimate in EUR, fara TVA si nu includ taxele de curierat.

Valoarea in LEI se va calcula la cursul BNR din ziua facturarii.

Produsele se vor ridica de la sediul nostru, sau se vor expedia in sistem ramburs prin curierat rapid, taxele de transport fiind suportate de destinatar.

Date contact :  
[andreescu@labor-soft.ro](mailto:andreescu@labor-soft.ro)  
[www.labor-soft.ro](http://www.labor-soft.ro)  
021 212 4875  
0741 038 830





Diametru 200 mm

Cod produs : 00 262 200

Pret : 60.35 EUR fara TVA (pretul nu include placa din portelan)

Bucati disponibile : 1

Produs marca Simax



Cu talpa si cioc, marcaj brun - linie scurta.

Cod produs : 2200190, volum 500 ml

Pret : 9.85 EUR fara TVA.

Bucati disponibile : 1

Produs marca Hirschmann



Din sticla DURAN, dop din plastic (PE), marcaj albastru

Cod produs	Volum	Pret EUR fara TVA	Bucati disponibile
2820190	500 ml	13,42	1

Produse marca Hirschmann



Clasa B, gradate pana la varf, marcaj brun, gradatie 0.01 mL, toleranta 0.015 mL, lungime 360 mm.

Cod produs : 1000140, volum 2 ml

Pret : 1.11 EUR fara TVA.

Bucati disponibile : 3

Produse marca Hirschmann



Clasa AS, gradate pana la varf, marcaj brun

Cod produs	Volum	Pret EUR fara TVA	Bucati disponibile
1100133	1	1,54	2
1100160	10 ml	1,99	2
1100165	20 ml	3,64	1
1100175	50 ml	5,67	3

Produse marca Hirschmann



Clasa AS, gradatie pana la varf, gradatie albastra, toleranta 0.01

Cod produs : 1101141, volum 2 ml

Pret : 1.40 EUR fara TVA.

Bucati disponibile : 1

Produs marca Hirschmann



Clasa AS. Nu sunt gradate pana la varf. (cu doua cote : gradatia de 0 sus si gradatia de 2 in apropierea varfului). Cu gradatii brune.

Cod produs : 1100241, volum 2 ml

Pret : 1.49 EUR fara TVA.

Bucati disponibile : 7

Cod produs : 1100260, volum 10 ml

Pret : 1.66 EUR fara TVA.

Bucati disponibile : 1

Produs marca Hirschmann





Clasa AS, cu o singura gradatie, bruna.

Cod produs	Volum pipeta	Pret/buc fara TVA	Bucati disponibile
<a href="#">1340129</a>	0,5 ml	1,74	2
<a href="#">1340133</a>	1 ml	1,52	2
<a href="#">1340155</a>	5 ml	2,44	2
<a href="#">1340160</a>	10 ml	2,76	2
<a href="#">1340170</a>	25 ml	3,17	1
<a href="#">1340173</a>	30 ml	4,54	4
<a href="#">1340175</a>	50 ml	4,26	1



Clasa AS, cu o singura gradatie, albastra.

Cod produs 1341133, volum 1 ml

Pret 1,38 EUR/buc fara TVA, bucati disponibile 6

Cod produs 1341140, volum 2 ml


Pret 1,89 EUR/buc fara TVA, bucati disponibile 5

Produse marca Hirschmann



Clasa AS, cu doua gradatii, brune.

Cod produs	Volum	Pret EUR, fara TVA	Bucati disponibile
<a href="#">1340233</a>	1 ml	<a href="#">1,98</a>	11
<a href="#">1340240</a>	2 ml	<a href="#">2,77</a>	10
<a href="#">1340260</a>	10 ml	<a href="#">3,61</a>	14
<a href="#">1340265</a>	20 ml	<a href="#">4,67</a>	2
<a href="#">1340275</a>	50 ml	<a href="#">5,73</a>	1



Clasa B, gradatie neagra, cu robinet din sticla, volum 50 ml, diviziuni 0.1 ml

Cod produs 3100175

Pret 21,34 EUR/buc fara TVA, bucati disponibile : 1

Produse marca Hirschmann



Clasa B, sticla DURAN®, gradatie neagra, cu robinet drept PTFE, volum 25 ml, diviziuni 0.05 ml, calibrata pt. EX, DIN ISO 385

[Cod produs 3100270](#)

[Pret 21,34 EUR/buc fara TVA](#), bucati disponibile : 4

Produse marca Hirschmann



Din sticla borosilicat 3.3, forma inalta, cu gradatii si cioc

Cod produs : 7213391, volum 600 ml

Pret : 2.58 EUR fara TVA, bucati disponibile : 15

Cod produs : 7210192, volum 800 ml

Pret : 6.51 EUR fara TVA, bucati disponibile : 7

Produse marca Hirschmann



Din sticla borosilicat 3.3, gat ingust cu margine bordurata

Cod produs : 7000186, volum 250 ml

Pret : 3.04 EUR fara TVA, bucati disponibile : 5


Cod produs : 7003387, volum 300 ml

Pret : 2.1 EUR fara TVA, bucati disponibile : 4

Cod produs : 7000193, volum 1000 ml

Pret : 7.46EUR fara TVA, bucati disponibile : 3

Produse marca Hirschmann



Din sticla DURAN borosilicat 3.3, gat ingust, cu slif NS 29/32

Cod produs : 7023385, volum 200 ml


Pret : 4.9 EUR fara TVA, bucati disponibile : 1

Cod produs : 7020186, volum 250 ml

Pret : 6.18 EUR fara TVA, bucati disponibile : 43

Produse marca Hirschmann





Din sticla DURAN, cu gat larg si margine bordurata, gradatie alba.  
Diametru gat 50 mm, inaltime 220 mm.

Cod produs : 7100193, volum 1000 ml

Pret : 7.09 EUR fara TVA, bucati disponibile : 3

Produse marca Hirschmann



Din sticla DURAN, prevazut cu tub pentru conectare la furtun

Cod produs : 2120144, volum 500 ml

Pret : 17.09 EUR fara TVA, bucati disponibile : 7

Cod produs : 2120154, volum 1000 ml

Pret : 31.80 EUR fara TVA, bucati disponibile : 3

Produse marca Schott



Din sticla DURAN. Vas de filtrare din sticla cu racord lateral vid din PP.

Cod produs : 212014445, volum 500 ml

Pret : 27.11 EUR fara TVA, bucati disponibile : 6

Produs marca Schott



Din sticla borosilicat 3.3, unghi de deschidere de 60°

Cod produs : 00 276 035

diametru 35 mm, capacitate 8 ml, lungime tija 35 mm, diametru tija 6 mm

Pret : 2.38 EUR fara TVA, bucati disponibile : 3

Cod produs : 00 276 075

Diametru 75 mm, capacitate 95 ml, lungime tija 75 mm, diametru tija 8 mm

Pret : 3.76 EUR fara TVA, bucati disponibile : 5



Din sticla borosilicat 3.3

Cod produs : 00 308 250

Volum 250 ml, inaltime 270 mm, fara slif

Pret : 9.04 EUR fara TVA, bucati disponibile : 4



Din sticla Duran, volum 250 mL, cu gat lung

Cod produs : 2171136

Pret : 4.8 EUR fara TVA, bucati disponibile : 4

Producator : Schott



Volum 500 mL, NS 29/32

Cod produs : 00 322 519

Pret : 3.71 EUR fara TVA, bucati disponibile : 3



Eprubete gradate, cu dop din PP, slif 14/23, gradatie albastra.  
Volum 25 ml. inaltime 215 mm

Cod produs : 6420170

Pret : 5.57 EUR fara TVA, bucati disponibile : 5

Producator : Hirschmann





Cod produs : 00 643 060

Margine slefuita, diametru 60 mm

Pret : 0.66 EUR fara TVA, bucati disponibile : 4

Cod produs : 00 643 150

Margine slefuita, diametru 150 mm

Pret : 0.76 EUR fara TVA, bucati disponibile : 8

# ANDREESCU LABOR & SOFT SRL

Stoc produse din sticla

Preturile sunt exprimate in EUR, fara TVA si nu includ taxele de curierat.

Valoarea in LEI se va calcula la cursul BNR din ziua facturarii.

Produsele se vor ridica de la sediul nostru, sau se vor expedia in sistem ramburs prin curierat rapid, taxele de transport fiind suportate de destinatar.

Date contact :  
[andreescu@labor-soft.ro](mailto:andreescu@labor-soft.ro)  
[www.labor-soft.ro](http://www.labor-soft.ro)  
021 212 4875  
0741 038 830

Adresa punct de lucru :

Str. Leonida, nr. 13, sector 2  
Bucuresti

